

## DEFIBRILLATOR ANALYZERS

### AED TESTER DELTA 1600

**NEW!**

- Unique, stand-alone instrument for testing all AEDs
- Multiple arrhythmia simulations
- Accurate energy measurements
- Automatic shock advisory testing
- Large backlit Graphics Display
- Li-Ion Rechargeable Battery



### DEFIBRILLATOR / TRANSCUTANEOUS PACEMAKER ANALYZER DELTA 3300

**NEW!**

**Free ProCom Print software  
of value \$200 included!**

- Precision instrument for testing defibrillators and external (Transcutaneous) Pacemakers
- Ideal for testing all defibrillators including AED's
- Waveform storage and Playback capability
- Perform Test, Save results and Print using ProCom Print software
- Large backlit LCD Display
- Compact & Lightweight (3lbs)
- Simple & user friendly operations.
- Li-Ion Rechargeable Battery
- Incredible value in comparison with competition products



**TEST.  
SAVE.  
PRINT.**

## EXTERNAL PACEMAKER ANALYZER

### EXPMT 2000

- Precision instrument for testing any External Pacemaker
- Menu driven and easy to operate.
- 12 Selectable Test Loads and 5 ECG Test Leads
- Performs Demand Sensitivity & Sensed and Paced Refractory tests.
- Large LCD Display
- Battery Operated



**MADE IN USA!**

## DIGITAL GAS FLOW METER

### DIGIFLOW 1000

- Precision instrument for gas flow measurements (Air, N<sub>2</sub>, O<sub>2</sub>, CO<sub>2</sub>, He)
- Large LCD Display
- Compact & Reliable
- Simple Menu Driven and User Friendly operation.
- Battery Operated



## INFUSION PUMP ANALYZER

### IPA 2000

- Single-channel accurate infusion pump flow tester with Netech's Patented Technology
- Measures Flow Rate, Volume, Occlusion
- Tests All Infusion Pumps, Including Syringe Pumps
- Extremely Accurate even at Low Flow Rates
- Large backlit Graphics LCD Display With Full Function Keypad
- Automatically Saves Test Results to Internal Flash Memory
- Uploads Stored Test Data to a PC Via communication port
- Menu Driven, Easy to Operate
- Mounts on Any IV Pole



## CONTACT US



**1-800-547-6557**



**Sales@NetechCorp.us**



**www.NetechCorp.us**



**@NetechBiomedical**



Serving Biomedical Industry  
Worldwide Since 1987

## INNOVATIVE BIOMEDICAL TEST INSTRUMENTS



**COMPACT. EASY TO USE.  
BEST VALUE.**



Contract Holder



REGISTERED



8a Certified

**WWW.NETECHCORP.US**

## ELECTRICAL SAFETY ANALYZERS

### SAFETY ANALYZER W/OUT PATIENT LEADS LKG 601

- Lightweight, Portable & Rugged design
- Large LCD display
- Accessible line fuse
- Calibrated test points
- 15 and 20 Amp capability
- Complies with AAMI or IEC 601-1 Test Loads




### SAFETY ANALYZER WITH PATIENT LEADS LKG 610

- Lightweight, Portable & Rugged design
- Large LCD display
- Point-to-Point testing
- Calibrated test points
- Patient Lead Isolation Test
- 15 and 20 Amp capability
- Complies with AAMI or IEC 601-1 test loads



### FULL-FUNCTION SAFETY ANALYZER MULTIPRO 2000

**NEW!**



- Quick and easy to operate
- Microprocessor controlled
- 12 Lead ECG and performance Waveforms simulation
- 20 Amp Rating
- Accurate resistance reading with Kelvin Cable
- Large Display for easy reading
- Communication Port
- Save test results and Print using ProCom software
- Complies with AAMI or IEC 601-1 test loads
- Free ProCom Software Included!

## PATIENT SIMULATORS

### ECG ARRHYTHMIA SIMULATOR MINISIM 1000 BASIC

- 12 lead ECG simulation with 14 user selectable waveforms.
- Menu driven operation
- Large LCD display
- Tactile Keypad
- Battery or AC operation



### FULL FUNCTION PATIENT SIMULATOR MINISIM 1000 ADVANCED

- 12 lead ECG simulation with 45 Arrhythmias
- ECG Performance waveforms
- Respiration / Apnea Simulation
- 2 channels IBP Simulation
- YSI 400 & YSI 700 Temperature probe simulation
- Menu driven operation
- LCD display
- Battery or AC operation



### EEG SIMULATOR MINISIM EEG

- Ideal simulator for sleep study & EEG monitor testing
- 5 channels patented cup electrode connectors
- Large LCD display
- Clean microvolt signals
- User Selectable amplitude & frequencies
- High accuracy
- Battery or AC operation



### CARDIAC OUTPUT SIMULATOR MICROSIM COS

- Baxter Edwards Catheter Compatible
- Ideal for simulation of: 4 blood flow rates, 2 blood temperatures, 2 injectate temperatures
- Power up self test
- Battery Operated



## PRESSURE METERS

### DIGIMANO 1000

- Portable and lightweight digital air pressure and vacuum meter
- Displays pressure in two different engineering units
- Ideal for calibration of IBP cuffs / sphygmomanometers
- Special lure-lock port for easy connection
- Battery operated



### DIGIMANO 2500

- Precision digital air pressure / vacuum meter
- Displays pressure in multiple engineering units
- Large backlit Graphics LCD Display
- Automatic Zero, and Peak Pressure Reading features
- Communication Port
- Real Time Clock
- Data Storage Capability
- Rugged Aluminum Case
- Rechargeable Li-ion battery



### BP-CALIBRATION KIT

- Accessory Kit for blood pressure calibration
- Can be used with Digimano 1000 or Digimano 2500 pressure meters
- Manual adjustable syringe pump for precise steady infusion of pressure
- Includes flexible silicone tubing with connector pack
- Customized hard carrying case



### UNIMANO 1000

- Precision meter for measurement of pressure and vacuum
- Designed for medical, laboratory and industrial applications
- Special lure-lock port for easy connection
- Displays pressure in multiple engineering units
- Fluid and Gas Compatible
- Auto Zero, Leak Detection & Peak Hold features
- Gauge and Differential pressure measurement
- Two Line, 16 Character LCD Display
- Battery Operated

